

MICROPOINT

qLabs® PT-INR Test Strips

REF QS-1-12 Contains: 12 test strips



REF Q-1

qLabs® ElectroMeter

REF QS-1-24 Contains: 24 test strips



REF Q-3

qLabs® ElectroMeter

For Self-Testing Use Only

INTENDED USE

The qLabs® PT-INR test system consists of an ElectroMeter and PT-INR Test Strips. The qLabs® test system is designed to provide quantitative measurements of Prothrombin Time (PT) and International Normalized Ratio (INR) in fresh capillary whole blood. Fresh venous blood or anticoagulated blood shall not be used.

The qLabs® PT-INR Test Strips are intended for in vitro diagnostic use by patients taking oral anticoagulants (blood thinners) who need to monitor the clotting time of their blood.

For self-testing use only.

INTRODUCTION

Prothrombin time (PT) is the test of choice for monitoring patients who are receiving oral warfarin therapy. The international normalized ratio (INR) is the recommended method for reporting PT results that are independent of PT methods. INR plays a critical role in maintaining the warfarin response within a therapeutic range such that it provides the efficacy of anticoagulation (blood thinning) while avoiding the risks of haemorrhage. The qLabs® PT-INR system can be used to monitor the INR levels of patients undergoing warfarin therapy.

TEST PRINCIPLE

The qLabs® PT-INR system uses a one-stage PT method to measure the patient's INR levels. After a drop of blood is added to the strip, the blood flows to the test zones where it reacts with reagents, initiating clot formation. As clotting proceeds, the qLabs® ElectroMeter detects the change of electric current across the clot, which is used to determine PT and INR results.

QUALITY CONTROL

The qLabs® system utilizes a number of internal quality methods to ensure proper operation. The built-in quality control of instrument automatically monitors critical conditions before and during the testing period. The onboard quality control of strip detects the signal characteristics of two channels. By identifying possible problems such as test strip defects and operational issues, to ensure the accuracy of the test results.

Failure of Quality Control test will result in the instrument displaying an error code. Please repeat the test utilizing a new test strip.

REAGENTS

Each test strip contains:

- Reagent channels: Recombinant human thromboplastin
- Both channels: Heparin neutralizing reagent

PRECAUTIONS & WARNINGS

- For in vitro diagnostic use only. Do not take internally.
- Follow proper infection control guidelines for handling all blood specimens and related items.
- Use fresh capillary blood ONLY.
- Never add blood to a test strip after the test has begun.
- Do not use strong repetitive pressure to collect the sample.
- Do not move the meter during a test.

The health status of the patient may affect the test. Please take this into consideration before making a therapeutic judgment based on the test results. Failure to do so may have serious consequences.

See the results section below for more information.

STORAGE & HANDLING

qLabs® PT-INR Test Strips can be stored at room temperature (below 32° C) or in the refrigerator at 2° C to 8° C until the expiration date. Do not freeze.

Store strips in their original foil pouch until ready to use.

If refrigerated, allow the sealed pouch to equilibrate to room temperature for 5 minutes before opening it for testing.

Use the test strip within 10 minutes of opening the foil pouch.

SAMPLE PREPARATION

Materials provided

- qLabs® PT-INR Test Strips
- CodeChip (for Q-3 instrument only)

Materials required (but not provided)

- qLabs® ElectroMeter
- Puncture-resistant container for medical sharps
- Alcohol Pads and Gauze
- 23-gauge or larger Lancet Device

⚠ Make sure the hand is warm. If not, warm the hand by washing in warm water or using a heating pad.

TEST PROCEDURE

⚠ When the meter powered on, refer to the User's Manual of qLabs® ElectroMeter to enter the Test Mode and prompt you to insert a test strip

1. **Insert a test strip into the test strip guide on the meter.** Remove a fresh test strip from its foil pouch. Insert it into the test strip guide so that the electrode end goes in first. On the light blue end of the strip you should be able to read the word "PT-INR" appearing from left to right.
2. **Input the Strip Code / strip Codechip number.**
 - 2.1 For Q-1 meter, input the **Strip Code** information labelled on the pouch. Then check the **Strip Code** to see if it is the same as the code on the pouch, correct once the code is wrong
 - 2.2 For Q-3 meter, input the strip **Codechip number** labelled on the pouch. Then insert the **Codechip** of the strip into the chip slot. qLabs will automatically check if it matches the strip **Codechip number** entered. If not, meter will display an error, and user needs to re-test by entering the right strip **Codechip number** or inserting the right **strip Codechip** to continue the

test.

⚠ Always match the Strip Code or Codechip number on the display with these on the strip pouch. Failure to do so may yield inaccurate results.


Strip Code




Codechip number



3. **Wait for the meter to warm up.** The ElectroMeter will warm up automatically for the test. When it is ready to perform a test, the ElectroMeter will beep and prompt user to apply a blood sample.
4. **Obtain a fingerstick blood sample.** It is important that you use the correct technique to obtain the right type and amount of blood sample. If the procedure is not followed, it can cause inaccurate results.



- 4.1 Increase blood circulation by:
 - Warming the hand with a heating pad or hand warmer
 - Gently massaging the finger
 - Holding the hand below the heart
- 4.2 Identify a site on the finger to puncture:
 - On one of the middle fingers of either hand
 - Near the top of the finger on either side
 - Away from any calluses or scars
- 4.3 Clean the selected area with 70% isopropyl alcohol, or an alcohol pad. Dry thoroughly with cotton or gauze.
- 4.4 Puncture the finger following the instructions for the lancet that you are using.
- 4.5 Apply gentle, continuous pressure until a large, hanging drop of blood (at least 10µL) forms. Do not use strong repetitive pressure to collect the sample.
- 4.6 Add the hanging drop of blood to the sample well of the test strip.



5. **Perform PT test.** After adding blood sample, the system will start test automatically. The test results will appear on the screen.
6. **Finish the test.** Discard the used lancet and test strip into a puncture resistant waste container. All blood samples should be regarded as potentially hazardous.

RESULTS

Since PT results are expected to vary with the test method, it is recommended that the same method must be used whenever doing routine patient monitoring. Each lot of test strips is calibrated to a reference lot that is traceable to the WHO International Reference Preparations. The International Sensitivity Index (ISI) for the system has been established as close to 1.

Normal Range:

Results for normal blood were determined by testing 120 subjects who were not taking anticoagulant medication. The ranges found were: INR: 0.70 -1.40.

Therapeutic Range:

Therapeutic ranges are determined for each patient individually by their clinical professional. While most recommendations are to be within an INR range of 2.00 to 4.50, values well below or well above that may be encountered.

Unexpected Results:

When the ElectroMeter displays a PT-INR result outside of the expected therapeutic range, it may or may not be due to an unusual clinical situation.

What causes unexpected results:

Hematocrit: The qLabs® system is validated to work reliably with blood having hematocrit values between 30% and 55%. Blood samples outside of this range may give unusual PT values.

Interfering antibodies: Conditions (such as Lupus) that produce antiphospholipid antibodies may interfere with the ability of blood to clot through the normal means.

Interfering metabolites: The qLabs® system is validated to work in the presence of unusually high concentrations of hemoglobin, bilirubin, or triglycerides (see Limitations Section below). Presence of these metabolites at concentrations above these limits may lead to long clot times.

Medications: Certain medications, including both prescription and over the counter, may interfere with oral anticoagulants, and may lead to an anomalous INR result.

Disease state: Certain medical conditions may interfere with anticoagulant therapy.

Diet: Oral anticoagulants may be sensitive to food, alcohol, and nutritional supplements.

What to do:


Whenever you encounter an unexpected result, please repeat the test with a fresh qLabs® test strip. If the result is seen a second time, please consult immediately with your health care professional and local distributor. Please don't take any decision of medical relevance without first consulting the appropriate health care professional. If any serious incident related to the instrument has occurred, please report it to us and the competent authority of the Member State in which you are established.

LIMITATIONS

- The qLabs® system is designed to use fresh capillary whole blood. Fresh venous blood or anticoagulated blood shall not be used.
- The qLabs® system is not affected by Heparin concentrations up to 1 anti-Xa units per mL of blood. This is true for both unfractionated heparin and low molecular weight Heparin.
- The drop of blood drop must be at least 10 µL in volume or a large hanging blood drop. Low sample volume will cause an error message.
- In vitro studies show no significant effect in blood samples containing up to 20 mg/dL of bilirubin, 500 mg/dL of hemoglobin (hemolysis), or 1500 mg/dL of triglycerides (lipemia).
- The qLabs® PT-INR Test Strips are validated to perform at temperatures in the range 10 to 35°C, and 10 to 90% RH (relative humidity). This includes a 10 minute out of pouch exposure of the strips at these conditions.

ACCURACY

Regression analysis of the qLabs® system Compared to the Werfen ACL-TOP 700
 $(y = 0.8052x + 0.436, n=101, r=0.922.)$



ADDITIONAL INFORMATION

If you have any questions regarding the use of this product, please contact the local distributor or Micropoint Technical Support by emailing customerservice@micropointbio.com or calling +86 755 21600849.

PERFORMANCE CHARACTERISTICS

Category	Performance Specification
Intended sample	Fresh fingerstick whole blood
Operating temperature range	10 - 35°C
Operating humidity range	10 - 90% RH
Out-of-pouch stability	10 minutes
Shelf life (room temp, in pouch with desiccant)	18 months
Measurable range	INR between 0.50 - 8.00
Accuracy	Reference to ACCURACY part
Precision	CV ≤ 5%
Factor II sensitivity	13%
Factor V sensitivity	48%
Factor VII sensitivity	45%
Factor X sensitivity	57%
Hematocrit range	30 - 55%
Time to results	30-100 seconds
Sample volume	≥ 10 µL

SYMBOLS EXPLANATION

Symbols	Explanation
	In vitro diagnostics
	Name and Address of Manufacturer
	European Authorized Representative
	CE Marking
	Temperature limitation
	Lot number
	Date of Manufacture
	Expiry Date
	Do not reuse
	Catalogue number
	Contains sufficient for n tests
	Caution! Read Carefully.
	Consult instructions for use

REVISION HISTORY

Rev. No	Description of Change	Revision Date
A1	New Release.	2020-06-30

MICROPOINT

Micropoint Biotechnologies Co., Ltd.
 3-5F, Building 1, Runheng Electronics Factory
 Liuxian 2 Road, Xinan Street, Baoan District
 518101 Shenzhen, China
 customerservice@micropointbio.com
 www.micropointbio.com
 Tel +86 755 21600849
 Fax +86 755 86673903

Obelis SA
 Bd. General Wahis, 53
 1030 Brussels, Belgium

 www.obelis.net
 Tel +32 2 732 59 54
 Fax +32 2 732 60 03

qLabs® and Micropoint® are registered trademarks of Micropoint Biotechnologies Co., Ltd.
 ©2020 Micropoint Biotechnologies Co., Ltd.
 All rights reserved. Printed in China.
 P/N 631-62002 Rev.A1 EN 2020-06-30