 MICROPOINT  M101-091001
mLabs® D-Dimer
Rapid quantitatiave microfluidic assay for the detection of D-Dimer

For Health Care Professional Use Only
--

INTRODUCTION

When fibrin in a blood clot is broken down by enzymes in a so called fibrinolysis, a series of Fibrin Degradation Products (FDP's) are formed. One of the unique FDP's is D-Dimer, which concentration provides information helpful for diagnosing thrombosis. D-Dimer test is of significant clinical use in making diagnostic decision for patients suspected of thrombotic disorders, such as deep venous thrombosis (DVT), pulmonary embolism (PE) or disseminated intravascular coagulation (DIC).

TEST PRINCIPLE

The mLabs® D-Dimer test is based on the immunoassay technology in a microfluic cartridge. The immunoassay technology relies on the inherent ability of an antibody to bind to the specific structure of a molecule. In the mLabs® D-Dimer microfluidic cartridge, two antibodies (reporter and capture) are carefully chosen so to have excellent specificity and sensitivity for D-Dimer. The reporter antibody with a fluorescent dye attached is loaded in the reaction zone of the cartridge, while the capture antibody is immobilized in the detection zone. As the sample flows through, the D-Dimer in the sample first encounters and binds to the reporter antibody. As the sample continues it's marching through the detection zone, the capture antibody captures the D-Dimer reporter complex by forming a sandwich structure and remains in the detection zone, while the reporter antibody without D-Dimer is washed away. The fluorescent signal in the detection zone is proportional to the D-Dimer concentration in the sample.

MATERIALS PROVIDED

- 25 Test cartridges
- Pipette tips
- 1 IC card
- 1 Package insert

SAMPLE COLLECTION & STORAGE

Always wear protective gloves and suitable lab coats when handling patient samples as they may potentially be infectious.

All samples should be regarded as potentially hazardous and/or contaminated.

- Citrated is the only coagulant used for mLabs D-Dimer test.
- Ensure that the collection tube is completely filled to maintain the correct anti-coagulant to blood ratio.
- Thoroughly mix the whole blood sample through gentle inversion (6-8 times) of the tube.
- Perform the whole blood test within 60 minutes after sample collection.
- Refrigerate whole blood samples after collection.
- Do not use syringes to collect test samples.

TEST SAMPLE STABILITY

Whole blood samples are stable at room temperature for a maximum of 60 minutes.

Whole blood samples, if processed and refrigerated at 2-8°C within 60 minutes after collection, are stable for 6 hours.

Severely hemolyzed samples should be avoided.

TEST KIT STORAGE

Refrigerate the D-Dimer test kit at 2-8 °C immediatly upon receipt.

Only remove the number of tests required from refrigeration.

The kit can be stored for ONLY ONE day after removing from refrigeration.

The D-Dimer test kit is stable at 2-8 °C until its expiration date.

The cartridge can be retumed to refrigeration for only two times.

TEST PROCEDURE

1. Sample Preparation

- The D-Dimer cartridge is only to be used together with immunometer.

- Equilibrate the pouched test cartridge to room temperature (20-30 minutes) before use.

- Thoroughly mix the collected whole blood sample by gently inverting (2-3 times) the tube before test.

2. Immunometer Preparation

- A single DataDrive (IC card) is provided along with each kit of cartridges.

- Insert the IC Card into the IC card port located in the right side of the meter.

- Press “DataDrive Installation” button on the main screen and select “Sample DataDrive”, then press “IC card”. The meter will automatically install DataDrive.

- Once the DataDrive has been uploaded into the Immunometer, the IC card can be removed for all subsequent related measurements.

- Store the IC card in a clean and dry location for future use.

3. Sample Addition

- Remove test cartridge from pouch and label it with the patient’s ID on the front with a permanent marker.

- Using the provided pipette, transfer 250 µL of sample into the inlet of the cartridge, dropwise.

- Do not place the pipette tip into the inlet during sample transfer as air bubbles may be generated. Two test mode can be choosed, respectively internal mode and external mode which are discribed clearly in meter’s user manual. For internal mode, place the cartridge into the cartridge holder after pressing holder eject button, then add fresh blood into the sample well and close the holder, the meter will automatically test the sample after a few minutes (less than 16min). For external mode, after removing the cartridge from the pouch, add fresh blood sample into the sample well and place the cartridge on a platform for 10-15minutes, then press the holder eject button to place the cartridge into the holder and close the holder for testing.

4. Performing Test and Reading Results

- Insert the test cartridge to the cartridge holder of the Immunometer.

- Press “Sample Test” from the main screen of the Immunometer.

- Select the assay and sample type if needed.

- Enter patient ID if needed.

- Press “Confirm” to start testing.

- The results will be displayed on the screen after test is complete.

Note: For different software and meter version test procedure may be different, user can refer to meter’s user manual sample test part to get more detailed operation information.

PERFORMANCE CHARACTERISTICS

- Analytical Sensitivity
D-Dimer = 50 ng/ml
- Measurable Range
D-Dimer = 50 to 10,000 ng/ml
- Cut-off
D-Dimer = 500 ng/ml

INTERFERENCE TESTING

Hemoglobin (up to 5 mg/mL), lipids (triolein up to 30 mg/mL), bilirubin (up to 0.15 mg/mL) added to citrated-plasma containing D-Dimer did not interfere with the recovery of D-Dimer. These substances also did not generate a positive response in the absence of D-Dimer. It is noted that severely hemolyzed specimens should be avoided. Hematocrit in between 30% and 55% has no significant effect on the recovery of D-Dimer.

EXPECTED VALUES

In a study carried out using apparently healthy individuals' samples, the 95th percentile values were below 500 ng/mL. It is recommended that each laboratory should establish its own reference range.

LIMITATIONS

Carefully inspect the mLabs® D-Dimer test pouch’s integrity before use. If the pouch is found to be tampered with, i.e. torn or punctured, do NOT proceed to use the test cartridge. Contact your local technical support immediately.

The mLabs® D-Dimer test kit is strictly for In Vitro Usage only. Instructions and procedures provided in this insert should be carefully adhered to.

The mLabs® D-Dimer test kit is not intended to be used as absolute evidence for PE and DVT. Obtained test results should be consulted with physician in addition with other test results.

All the provided items in the test kit are for single usage application and should be properly discarded after usage as inserted test samples may potentially be infectious.

QUALITY CONTROL

The mLabs system has a number of on-board quality control functions integrated into the Meter and the cartridge and therefore there is no need to run quality control tests with quality controls. However, Micropoint has provided an optional quality control for the mLabs system. This control is provided to help meet the regulatory requirements applicable to your facility. Please contact your local distributor or call Micropoint Bioscience customer service at +1 408-588-1682 or email to customerservice@micropointbio.com.

REFERENCES

S.Z. Goldhaber. Pulmonary embolism. New England Journal of Medicine. 1998. 339: pg 93-104.

P.S. Wells, D.R. Anderson, M. Rogers et al. Evaluation of D-Dimer in the diagnosis of suspected deep-vein thrombosis. New England Journal of Medicine. 2003. 349. pgs 1227-1235.

P.F. Fedullo and V.F. Tapson. The evaluation of suspected pulmonary embolism. New England Journal of Medicine. 2003. 349. pgs 1247- 1256.

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©2019 Micropoint Bioscience, Inc. All rights reserved.Specifications are subject to change without notice. mLabs® and Micropoint® are registered trademarks of Micropoint Bioscience, Inc. Printed in China.